

# Detection of PARP-induced ADP-ribosylation using the CLARIOstar® microplate reader

- PARP enzyme activity determined using a chemiluminescence assay
- The CLARIOstar® luminescence readout is reproducible and linear over a wide enzyme concentration range
- The kinetic parameter  $K_M$  and the  $IC_{50}$  value for an enzyme inhibitor were calculated

## Introduction

PARP (Poly(ADP-ribose) polymerase) family enzymes are involved in the regulation of transcription, DNA repair, and chromatin remodeling.<sup>1</sup> These enzymes use nicotinamide adenine dinucleotide (NAD) as a substrate to build poly(ADP-ribose). Due to various links to diseases,<sup>2,3</sup> PARP enzymes are targets for pharmaceutical drug development.

In this application note we describe the use of a chemiluminescent assay to determine PARP activity on the CLARIOstar multimode microplate reader. The assay allows kinetic analysis of PARP enzymes and evaluation of inhibitor potency.

## Assay Principle

PARP activity is followed *in vitro* by detecting the incorporation of biotinylated ADP-ribose as a consequence of either enzyme target protein modification or auto-modification. The reaction principle is shown in Fig. 1.




Fig. 1: PARP chemiluminescent assay principle

Hexahistidine-tagged PARP enzyme or protein substrate is immobilized on Ni<sup>2+</sup>-chelating microplates. The reaction is started by adding biotinylated NAD<sup>+</sup>. The PARP enzyme uses the NAD<sup>+</sup> to synthesize biotinylated poly(ADP-ribose). This polymer is either added to the PARP enzyme itself or transferred to a protein substrate on the microplate (histone). After a washing step streptavidin-conjugated horseradish peroxidase is added and will bind to the biotinylated poly(ADP-ribose).<sup>4</sup>

After adding a substrate to the horseradish peroxidase, chemiluminescence is released and can be measured in the CLARIOstar.

## Materials and Methods

- Ni-NTA-coated, opaque, white 96-well microplates (5-PRIME, 2400730)
- Biotinylated NAD<sup>+</sup> (Trevigen, 4670-500-01)
- Streptavidin-conjugated horseradish peroxidase (Jackson ImmunoResearch, 016-030-084)
- SuperSignal West Pico chemiluminescent substrate (Thermo Fisher Scientific, 11513450)
- CLARIOstar microplate reader from BMG LABTECH (Fig. 2)

All standard chemicals and disposables were obtained through normal distribution channels.



Fig. 2: CLARIOstar multimode microplate reader from BMGLABTECH

## Enzymatic reactions

Hexahistidine-tagged PARP enzyme or protein substrate was immobilized on Ni<sup>2+</sup>-chelating plates. ADP-ribosylation reactions were started by addition of NAD<sup>+</sup> (2 % biotinylated) at 20°C. Reactions were stopped by addition of 7 M guanidine hydrochloride. Plate wells were washed with reaction buffer, incubated for 30 minutes with TRIS-buffered saline containing 0.02 % Tween-20 (TBST) and 1 % (w/v) BSA, and washed with TBST. After incubation with streptavidin-conjugated horseradish peroxidase (0.5 µg/ml) another washing step was done. After adding SuperSignal West (50 + 50 µl, undiluted) as substrate for the peroxidase chemiluminescence was detected in the CLARIOstar microplate reader using the following instrument settings.


## CLARIOstar instrument settings

All measurements (linear range check,  $K_M$  determination and inhibitor dose-response) were done in endpoint mode.

Optic: top optic used  
 Measurement interval time [s]: 1.00  
 Presetname: Enliten ATP  
 Emission: full range (no filter)  
 Gain: needs to be adjusted prior the measurement  
 Focal height: needs to be adjusted prior the measurement


## Results and Discussion

For validation of the ADP-ribosyltransferase assay the linear range of signals obtained by a dilution series of biotin-ADP-ribosylated enzyme was determined (Fig. 3). The results show that the signal is linear over a wide range of ADP-ribosyl concentrations.




**Fig. 3:** A dilution series of a PARP-family enzyme under assay conditions, illustrating the linear range of the signal. The insert zooms into the low nM concentration range.

The kinetic parameters of a PARP enzyme family member were determined using initial reaction rates. Independent experiments showed that the biotin moiety linked to the co-substrate had no influence on the reaction kinetics (results not shown).



**Fig. 4:** Rates plot of the NAD<sup>+</sup>-dependent ADP-ribosylation catalysed by a PARP-family enzyme.

Knowledge of  $K_M$  allowed the determination of inhibitor dose-response curves and experimental parameters ( $IC_{50}$ ).



**Fig. 5:** Dose-response curve for inhibition of a PARP-family enzyme with a clinical PARP inhibitor (Olaparib).

## Conclusion

An ADP-ribosylation assay of PARP enzymes carried out in a CLARIOstar microplate reader shows signal linearity over a wide range of enzyme concentrations (0.015 to 300 nM). The assay allows enzyme characterization and calculation of different parameters that are important for the development of drug-like enzyme inhibitors.<sup>5,6</sup>

## References

- Hottiger, M.O., Hassa, P.O., Lüscher, B., Schüler, H., and Koch-Nolte, F. (2010). *Toward a unified nomenclature for mammalian ADP-ribosyltransferases*. Trends Biochem Sci **35**, 208-219.
- Riffell, J.L., Lord, C.J., and Ashworth, A. (2012). *Tankyrase-targeted therapeutics: expanding opportunities in the PARP family*. Nature reviews Drug discovery **11**, 923-936.
- Curtin, N.J., and Szabo, C. (2013). *Therapeutic applications of PARP inhibitors: anticancer therapy and beyond*. Molecular aspects of medicine **34**, 1217-1256.
- Langelier, M.F., Planck, J.L., Servent, K.M. and Pascal, J.M. (2011). *Purification of human PARP-1 and PARP-1 domains from Escherichia coli for structural and biochemical analysis*. Methods Mol Biol **780**, 209-26.
- Franzini et al. (2015) *Identification of structure-activity relationships from screening a structurally compact DNA-encoded chemical library*. Angew Chem Int Ed Engl **54**(13): 3927-3931.
- Ekblad et al. (2015) *Towards small molecule inhibitors of mono-ADP-ribosyltransferases*. Eur J Med Chem **95**: 546-551.

Germany: **BMG LABTECH GmbH** Tel: **+49 781 96968-0**

Australia: **BMG LABTECH Pty. Ltd.** Tel: **+61 3 59734744**

France: **BMG LABTECH SARL** Tel: **+33 1 48 86 20 20**

Japan: **BMG LABTECH JAPAN Ltd.** Tel: **+81 48 647 7217**

UK: **BMG LABTECH Ltd.** Tel: **+44 1296 336650**

USA: **BMG LABTECH Inc.** Tel: **+1 877 264 5227**

Internet: **www.bmglabtech.com** **applications@bmglabtech.com**