

HTS Instrument Discovers Low Affinity Inhibitors of the Inositol Phosphate (IP) Signaling Pathway

EJ Dell¹, JL Tardieu², and F Degorce²,

¹ BMG LABTECH GmbH, 77799 Ortenberg, Germany

² Cisbio Bioassays, 30200 Codolet, France

Application Note 222, Rev. 02/2012

- Several low affinity drug 'hits' were only found by the PMT-based reader and not the CCD-based reader
- PHERAstar FS performs more than twice as fast as the CCD-based reader
- PHERAstar FS shows superior assay quality parameters such as Z', DeltaF% and assay window

Introduction

G-protein coupled receptors (GPCRs) are transmembrane proteins which play a key role in the signal transduction of extracellular stimuli. GPCRs are associated with a complex assembly of intracellular proteins regulating a large variety of downstream effectors. The production of inositol 1,4,5-triphosphate (IP₃) is one such second messenger, which is produced in response to the activation of G_q-coupled receptors. However, IP₃'s very short half-life make its assessment too challenging for drug screening assays and the monitoring of calcium release, triggered by IP₃, has been extensively used as a downstream readout of this signaling pathway.

An alternate way to monitor IP₃ is to measure the accumulation of inositol monophosphate (IP₁), which is a downstream metabolite of IP₃.




Fig. 1: Biosynthesis of IP₁.

Monoclonal antibodies raised against IP₁ lead to the optimization of a homogeneous time-resolved fluorescence assay (HTRF[®]), taking advantage on the fact that IP₁ is stable and accumulates in cells. The HTRF[®] IP-One assay has been compared to existing methods and has been shown to lead to similar compound potency data. Moreover, the end point accumulation of IP₁ allows for the discrimination of slow acting compounds that remain unseen by calcium sensing. The HTRF[®] IP-One assay also

allows for the characterization of inverse agonists by the quantification of constitutively active GPCRs, which is impossible via measurement of calcium release. Lastly, HTRF[®] IP-One detection confers superior assay robustness and much lower false positive rates compared to calcium detection. The IP-One assay developed by Cisbio Bioassays uses their proprietary HTRF[®] technology (Figure 2).

Assay Principle




Fig. 2: IP-One HTRF[®] Assay Principle

The assay uses a monoclonal antibody that specifically recognizes IP₁ and it is based on a competition format in which the intracellular accumulation of IP₁ inhibits the fluorescence resonance energy transfer (FRET) signal between the HTRF[®] donor and acceptor. An IP₁ calibration curve can estimate the IP₁ concentration accumulated in cells as a function of the compound concentration.

Materials and Methods

- Next generation PMT based HTS microplate reader, PHERAstar FS
- HTS CCD-based microplate imager from a different vendor




Fig. 3: 1536-well microplate layout

Black Greiner 1536-well microplates were used with an assay volume of 5 μ L. As shown in Figure 3, 1408 different compounds were pipetted into the first 44 columns of a microplate, with positive (POS, n=64) and negative (NEG, n=64) controls in the last 4 columns.

The PHERAstar FS is equipped with a high power pulsed nitrogen laser emitting at 337 nm, as well as a dedicated Simultaneous Dual Emission (SDE) direct photon counting time-resolved fluorescence mode. When exciting the terbium (Tb) donor molecule, the laser is superior to a broadband xenon flash lamp. The laser's energy emission takes advantage of the higher molecular extinction coefficient of the terbium cryptate peaking around 337 nm, compared to europium (Eu).

Results and Discussion

Results were evaluated by the HTRF[®] ratio of the two emission wavelengths (Em 665 nm / Em 620 nm) and hits are shown as peaks in a surface graph representing the entire 1536 well microplate.




Fig. 4: HTRF[®] ratios obtained for the IP-One assay with the PHERAstar FS




Fig. 5: HTRF[®] ratios obtained for the IP-One assay with a CCD-based HTS Device

Table 1: Speed and assay quality are compared between the PHERAstar FS and the CCD camera based HTS reader.

	PHERAstar FS	CCD based HTS reader
Read Times (1536)	53 sec	2:12 min
Assay Window	6:1	2:1
Delta F%	490	76
Z' Value	0.70	0.24

Conclusion

The IP-One HTRF[®] assay from Cisbio was performed on two different HTS microplate readers with different detection technologies. Low affinity compounds which were not discovered with a leading HTS CCD camera based imaging microplate reader were readily resolved with the PMT based PHERAstar FS.

This next-generation HTS reader, the PHERAstar FS from BMG LABTECH, represents a new choice for HTS screening assays.



Fig. 6: BMG LABTECH's multidetection microplate reader PHERAstar FS

Germany: **BMG LABTECH GmbH** Tel: **+49 781 96968-0**
 Australia: BMG LABTECH Pty. Ltd. Tel: +61 3 59734744
 France: BMG LABTECH SARL Tel: +33 1 48 86 20 20
 Japan: BMG LABTECH JAPAN Ltd. Tel: +81 48 647 7217
 UK: BMG LABTECH Ltd. Tel: +44 1296 336650
 USA: BMG LABTECH Inc. Tel: +1 877 264 5227
 Internet: **www.bmg-labtech.com** **applications@bmg-labtech.com**