

Protein-Ligand binding measurements using fluorescence polarization

Valentina Speranzini¹, Alexander Fish² and Andrea Mattevi¹

¹Lab of Structural Biology, University of Pavia, Dept. Biology and Biotechnology "L. Spallanzani", 27100 Pavia, Italy

²Division of Biochemistry and Center for Biomedical Genetics, Netherlands Cancer Institute, Amsterdam, The Netherlands

Application Note 256, Rev. 10/2014


- Fluorescence polarization assay used to probe protein-ligand interactions
- Binding affinities are consistent with activity assays
- Technique sensitivity allows analysis of a wide range of concentrations

Introduction


Human Lys-specific Histone Demethylase 1 (LSD1) promotes the demethylation of mono- or di-methyl-Lys4 on histone H3¹: as chromatin remodelling enzyme, LSD1 is found in complex with several partners, one of which being the REST Co-repressor 1 (CoREST1)². We used the fluorescence polarization (FP) as a biophysical tool to analyse the binding properties of LSD1-CoREST1 (LC1) hetero-dimer to H3-derived peptides (Fig. 1), that have been conjugated at the C-terminal with the fluorescent label TAMRA (5-Carboxytetramethylrhodamine).

Assay Principle

Direct Binding




(A) Rapid Tumbling
Low Polarization




(B) Slow Tumbling
High Polarization

Competition Assay



(C) Most complex is fluorescent
Highest polarization



(D) Most complex is non-fluorescent
Lowest polarization

Fig. 1: Assay principle for monitoring protein-ligand interaction using fluorescence polarization

(A) H3 labeled peptides in solution show a rapid rotation that causes dispersion in different direction of signal of the incident polarized light: the resulting emitted polarization signal is low.

(B) Upon binding to LC1, the high-molecular weight complex results in a slow rotation of the fluorescent molecule leading to high polarization values. (C,D) Non-tagged peptides, in orange, compete for LC1 active site. As competitor concentration increases, more freely labelled peptide remains in solution. This results in a low FP value, comparable to results in (A).

Materials and Methods

- black, 384-well, low flange microplates from Corning, UK
- TAMRA-conjugated peptides provided by NKI Protein Facility (Amsterdam)
- LC1 protein was recombinantly produced, according to reference 1
- CLARIOstar® multimode microplate reader from BMG LABTECH, Germany (Fig. 2)



Fig. 2: CLARIOstar multimode microplate reader from BMG LABTECH

Experimental setup

1 nM labelled peptides were mixed to assay buffer: 15 mM KH₂PO₄, pH 7.2, 5% glycerol, 1 mg/ml BSA. Mix was then distributed equally across the plate. For each experiment, triplicates were prepared: starting concentration of protein LC1 was 4 μM, and then this was serially diluted 1:1 using the mix in the other wells. For competitive assays, constant enzyme concentration (around K_d) were mixed with 1 nM labelled ligand in the same assay buffer as above. Competitors were then titrated in using serial 1:1 dilutions as described for direct binding assay. The typical concentration range started from 10 μM.

CLARIOstar instrument settings

Detection Mode: Fluorescence Polarization, Endpoint Mode
No. of flashes: 50
Filters: Ex: 540-20/LP565/ Em: 590-20
Target temperature: 25°C

Focus and Gain adjustments for both channels were set on a reference well containing just labelled ligand. It was chosen to set the adjusted polarization to a reference value of 35 mP. Typical adjustment values for gain were between 1600 and 2000.

Results and Discussion

We first devised a direct binding experiment using three different H3-derived peptides, all previously described as substrates and/or inhibitors of LC1 complex^{3,4}. On a relative scale, results obtained are consistent with activity measured before by the working group. The difference in the peptides consisted in different lengths and/or different amino-acid compositions as listed below:

- meH3 (pep01) = H3 N-terminal tail, methylated on Lys4, 21 residues
- H3K4M (pep02) = H3 N-terminal tail, bearing the mutate residue in position 4, 21 residues
- SN12 (pep03) = Transcription factor SNAIL, an H3 analogue, consisting of 12 residues (N-terminal)

As shown in figure 3, all three peptides do bind to LC1.




Fig. 3: Direct binding assay for LC1 and H3-derived peptides. In blue (pep01) H3 tail, 21 aa, methyl-Lys. In orange (pep02) H3 tail, 21 aa, K4M. In purple (pep03) SNAIL short peptide, 12 aa. Different amplitudes in the curves might be due to different protein-ligand conformations in solution.

Acquired data was analysed using Mars Software and graphically edited in GraphPad. Association constants were calculated modifying the equation for fluorescent anisotropy (see equation 1) using constant label concentration as constrain.^{5,6}

$$B = \frac{L_T + K_d + R_T - \sqrt{[(L_T + K_d + R_T)^2 - 4L_T R_T]}}{2}$$

Eq. 1: Fitting for Fluorescence Anisotropy. B= binding, Lt = total ligand concentration, Rt= total "receptor" concentration

To further assay the binding of histone tail derivate peptides and to confirm association parameters, we proceeded using a competition approach: for this set of experiments a new test mix was prepared, as described above. We titrated in different potential competitors, using the same enzyme concentration, which depended on Kd values obtained from direct binding assays. As shown below (Fig. 4), we were able to assess affinity of different potential binders using a different approach but still in fluorescence polarization.




Fig. 4: Competitive Assays for peptide binding to LC1 complex. In blue the reference direct binding curve is shown. In red and green are potential binders to the active site of LC1 that show a similar affinity for the enzyme.

Table 1 enlists all molecules used in both direct and competitive assays with the resulting affinities values (Kd).

Ligand	Kd	Type of Assay
Pep01	107.3 nM ± 7.6 nM	Direct Binding
Pep02	76.4 nM ± 10.4 nM	Direct Binding
Pep03	81.9 nM ± 16.7 nM	Direct Binding
Pep A	147.1 nM ± 35.4 nM	Competition
Pep B	316.6 nM ± 59.9 nM	Competition

Tab. 1: Summary of association constants for LC1-peptides binding. For pep01-02-03, the different constants are to be considered coherent, as differences might be due to the sensitivity of the measurement and in the mix preparation.

Conclusion

Fluorescence Polarization has proven to be a reliable technique to probe protein-ligand interaction in real-time and in solution. Both, direct and competitive assays, are useful tools to analyse binding properties. Next to the CLARIOstar very similar results were also obtained on the PHERAstar FS microplate reader.⁷

References

- Forneris et al. (2007). *JBC* 282(28):20070-4
- Yang et al. (2006). *Mol Cell* 23:377-87
- Karytinis et al. (2009). *JBC* 284(26):17775-82
- Lin et al. (2010). *EMBO J* 29:1803-16
- Jameson & Mocz (2005). *Methods Mol Biol* 305:301-22
- Mocz (2006). *J Fluoresc* 16:511-524
- Speranzini et al. (2015). *PNAS* 112:2752-7

Germany:	BMG LABTECH GmbH	Tel: +49 781 96968-0
Australia:	BMG LABTECH Pty. Ltd.	Tel: +61 3 59734744
France:	BMG LABTECH SARL	Tel: +33 1 48 86 20 20
Japan:	BMG LABTECH JAPAN Ltd.	Tel: +81 48 647 7217
UK:	BMG LABTECH Ltd.	Tel: +44 1296 336650
USA:	BMG LABTECH Inc.	Tel: +1 877 264 5227
Internet:	www.bmglabtech.com	applications@bmglabtech.com