

# Using CLARIOstar® fluorescence polarization detection to perform hERG Predictor™ assay

Carl Peters, BMG LABTECH, Cary, NC

- Predictor™ hERG Fluorescence Polarization Assay Kits represent an important test of cardiotoxicity
- Using the CLARIOstar® microplate reader you can expect large  $\Delta mP$  and Z' values, indicative of robust assay performance


## Introduction

Acquired long QT syndrome is a drug induced form of cardiac arrhythmia that is characterized by fainting and even sudden death as a result of ventricular fibrillation. The majority of drugs associated with this syndrome cause their effects by interacting with the potassium ion channel Kv11.1, a product of the human Ether-á-go-go-Related-Gene (hERG). Due to the cardiac associated issues it has been recommended that drug candidate interactions with hERG be evaluated in a preclinical environment.

The Predictor™ hERG Fluorescence Polarization Assay Kit was designed to fulfil the need to screen small molecules for the interaction with hERG and thus potential for cardiotoxicity<sup>1</sup>. Since it is based on fluorescence polarization (FP) it can detect hERG interactions without the need for radioactively labelled ligands while exhibiting a high correlation with patch clamp based techniques. Here, we demonstrate the performance of this assay using the CLARIOstar® microplate reader to assess the binding capability of several test compounds on a single microplate.

## Assay Principle

The Predictor™ hERG FP Assay uses a tracer which is a fluorescently labelled hERG channel ligand (Fig. 1).



**Fig. 1: Predictor hERG Assay Principle.** In the absence of a competing compound a high FP value will be observed. Competitors displace the tracer leading to a low fp value.

When the tracer is associated with the membrane bound hERG channel protein the rotation is slowed thus producing a high FP value. A compound that binds to the hERG channel protein will displace the tracer. When the tracer is free in solution it tumbles rapidly and thus exhibits a low FP value.

## Materials and Methods

- Predictor™ hERG Assay Kit (Life Technologies – PV5365)
- Test compounds (Sigma)
- Black, 384-well, low volume microplates (Corning – 3677)
- CLARIOstar multimode microplate reader from BMG LABTECH



**Fig. 2:** CLARIOstar multimode microplate reader from BMG LABTECH

Kit reagents were thawed and prepared according to the assay instructions and test compounds were dissolved in 100% DMSO. Dilution series for reference and test compounds were prepared using 15 point and 13 point series respectively with 3 fold dilutions. Dilutions were initially prepared at 100X with reference compound dilutions in Assay Buffer and test compound dilutions in 100% DMSO. A 4X intermediate dilution was subsequently prepared in Assay Buffer.


An assay plate was prepared by mixing the appropriate reagents, as indicated in kit instructions, to prepare 20  $\mu$ L each of Buffer Blank, Assay Blank, Free Tracer Control, Negative Control, Positive Control, E-4031 Titration (reference compound) and Test Compounds Titrations. The plate was then covered, protected from light and incubated for 2 hours at 20-25°C.

### Instrument settings

Measurement Method: Fluorescence Polarization, Endpoint Mode  
Filter Settings: 540-20 / LP 565 / 590-20  
Settling Time: 0.1  
Number of flashes: 200  
Focus and gain: adjusted prior to measurement  
Target mP: set to 50 mP for free tracer


## Results and Discussion

A spectral scan was performed on the Predictor™ hERG Tracer Red (Figure 3). This confirmed that the appropriate filters had been selected for detection of this assay.



**Fig. 3:** Excitation and Emission Spectra of Predictor hERG Tracer Red.


The assay performance could be assessed by comparing the positive and negative controls (Fig. 4).



**Fig. 4:** FP values for positive control (PC) and negative control (NC) of the hERG assay.


The positive control contains 30  $\mu\text{M}$  E-4031 along with tracer and hERG membranes while the negative control does not contain a hERG inhibitor. Comparing these controls show that performing this assay on the CLARIOstar® provides a very good assay window ( $\Delta\text{mP} > 130$ ) while a  $Z'$  value of 0.752 was achieved.

Assay performance is further exemplified by the reference compound titration curve (Figure 5).



**Fig. 5:** Reference compound titration curve.

Figure 6 shows the results from 2 test compounds one, HC-067047, was previously reported to interact with the hERG channel ( $\text{IC}_{50} = 370 \text{ nM}$ ) while the other, ML218, was reported to be a potent T-Type calcium channel inhibitor<sup>3</sup>.



**Fig. 6:** HC-067047 and ML218 titration curves. Test compound titration curves were prepared with a starting concentration of 11  $\mu\text{M}$  and employed in the hERG assay. The results from HC-067047 treatment conformed to a 4-parameter fit-curve ( $R^2 = 0.996$ ) with an  $\text{IC}_{50}$  of 68.4 nM. The data from ML218 did not conform as well to a 4-parameter fit-curve ( $R^2 = 0.948$ ) and predicted a high  $\text{IC}_{50}$ .

## Conclusion

The CLARIOstar® exhibits robust detection of the Predictor™ hERG assay as exemplified by the  $Z'$  values and  $\Delta\text{mP}$  calculated from controls. Furthermore, the reference compound and the test compounds performed as would be predicted based on previous results.

## References

- 1) Piper, D.R., *et al.* (2008) Development of the Predictor hERG Fluorescence Polarization Assay Using a Membrane Protein Enrichment Approach. *Assay Drug Dev. Technol.* **6**, 213-223
- 2) Everaerts, W., *et al.* (2010) Inhibition of the cation channel TRPV4 improves bladder function in mice and rats with cyclophosphamide-induced cystitis. *Proc. Natl. Acad. Sci. USA* **107**, 19084-19089
- 3) Xiang, Z., *et al.* (2011) The Discovery and Characterization of ML218: A Novel, Centrally Active T-Type Calcium Channel Inhibitor with Robust Effects in STN Neurons and in a Rodent Model of Parkinson's Disease. *ACS Chem. Neurosci.* **2**, 730-742

Germany: BMG LABTECH GmbH Tel: +49 781 96968-0

Australia: BMG LABTECH Pty. Ltd. Tel: +61 3 59734744

France: BMG LABTECH SARL Tel: +33 1 48 86 20 20

Japan: BMG LABTECH JAPAN Ltd. Tel: +81 48 647 7217

UK: BMG LABTECH Ltd. Tel: +44 1296 336650

USA: BMG LABTECH Inc. Tel: +1 877 264 5227

Internet: [www.bmglabtech.com](http://www.bmglabtech.com) [applications@bmglabtech.com](mailto:applications@bmglabtech.com)